



VACCINE FRIDGE TEMPERATURE CHART

Fridge: _____
Month: _____
Year: _____

Check and record the daily
min & max vaccine fridge
temperatures with an 'X',
then reset the thermometer.

www.skoep.com.au

°C \ Date	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
>15																																
15																																
14																																
13																																
12																																
11																																
10																																
9																																
8																																
7																																
6																																
5																																
4																																
3																																
2																																
1																																
0																																
<0																																

The vaccine fridge should operate between 4°C and 8°C. Remember – STRIVE FOR FIVE

If the vaccine fridge is outside of this temperature range, please isolate the vaccines and contact your State/Territory health authority for advice on the National Immunisation Program vaccines and contact the manufacturer/supplier for privately purchased vaccines.



Designed by SKOEP, suppliers of vaccine fridges, data loggers and min/max thermometers

www.skoep.com.au